

# XF-7800/7700/7600/7500

## Color Doppler Ultrasound Diagnosis System Series

Xianfeng 7 Series Color Doppler system has a unique shape of elegant atmosphere and the latest generation of application digital color Doppler diagnostic system which can be widely applicable to diagnose abdomen, obstetrics, gynecology, cardiac, mammary gland, musculoskeletal and peripheral vessel.



Optional probes

3.5MHz Convex



7.5MHz Linear



6.5MHz Trans-vaginal



5MHz Micro-convex



Volume Probe



### Parameters:

- 15" HD LED Displayer. Ergonomic design, All-round display angle adjustment
- Ultra-wideband imaging/THI/Trapezoidal imaging/ Exact SAR adaptive speckle noise suppression techniques.
- CFM mode / PDI mode / PW , Real time triplex Imaging.
- 6 different languages System.
- Rich measurement package which can satisfy clinical needs.
- support convex, linear, trans-vaginal, micro-convex, phased array and 4D probes.
- 3D/4D Optinal function.
- One-click optimization function.
- Image and video one-click fast storage.

